

Curtis L. Barmby, D.D.S.

Diplomate, American Board of Prosthodontics

Gail E. Frick, D.M.D.

Member of American College of Prosthodontists

HOW TO GET THE MOST OUT OF YOUR NEW MUCOSOFT™ RELINE

Your dentist has treated you with one of dentistry's advanced materials designed to improve the comfort and stability of your denture. Like so many other things in life, if you pay attention to proper maintenance, you'll enjoy this treatment over a long period of time. **So here are a few HELPFUL HINTS on how to care for your denture.**

DO place a washcloth on the bottom of your sink whenever washing your dentures so that, if you accidentally drop the denture, its fall will be "cushioned".

DON'T scrub the MUCOSOFT™ soft cushion liner with any kind of brush. Instead use cold running water and gently wipe it with cotton or soft gauze, making sure to remove all food and debris. **Only the hard plastic sections of the denture can be cleaned with a soft bristle toothbrush and a paste denture cleaner or toothpaste.**

DO soak your denture on a daily basis, preferably overnight in a denture solution like Polident®, Efferdent®, or similar product. **NEVER** use strong bleach like Clorox®, vinegar, etc. If you do, the color of the MUCOSOFT™ liner will fade, and its useful life may be shortened

ALWAYS keep your denture in water if it is out of the mouth for an extended period, like overnight. **DO** have your denture thoroughly examined and professionally cleaned by your dentist **every 3 months.**

DO inspect your denture **every day.** If the soft liner shows signs of cracking, breaking, loosening or is damaged in any way, stop wearing it and see your dentist so the problem can be fixed.

STOP wearing your denture if there are any signs of irritation or gum soreness. See your dentist as soon as possible to determine the cause.